
Gebruiksaanwijzing - LED Gezichtsmasker
User Manual - LED Face Mask
Manual d'utilisation - Masque LED pour le visage
Gebrauchsanweisung - LED Gesichtsmaske

Introductie

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Introduction

Thank you for choosing our product! Always follow the safety instructions for safe use. If there are any questions about the product please feel free to contact us via email: info@good-4you.nl.

Intended Use

- This product is intended only for the purposes described in the manual.
- The supplier is released from liability for any injury or damage caused by using the product beyond its intended use.
- Always observe the following precautions when using 'the LED Face Mask', hereinafter referred to as 'the product'.

Warnings and Safety Instructions

Make sure that the packaging material does not get into the hands of children. They could choke on it!
Read the instructions for use carefully before use and keep them in a safe place!

General

- If this is your first time using the product, it is important to read the manual and instructions carefully. Note the risks and consequences associated with misuse of the product. Misuse of the product may result in personal injury and damage to the product or products connected to it.
- Use the product only for the purposes indicated in this manual.
- If the product is malfunctioning, dropped or damaged, do not use the product and consult Customer Service.
- The product is designed and intended for private use. Do not use the product for commercial purposes.
- Do not leave the product unattended when in use. When not in use, turn the product off

Persons

- This product is not suitable for persons (including children) with physical or mental disabilities, or for persons who lack experience or knowledge of the product. Use of the product is only permitted under supervision by a person responsible for the user's safety. Use and store the product out of the reach of children.
- We recommend using this product only if you are 16 years of age or older.
- Avoid using this product on the same day as products containing light-sensitive ingredients.
- Do not use this product in combination with implanted electronic devices, such as artificial hearts, artificial lungs, etc. If in doubt, consult a physician first.
- Do not use this product if you have a genetic eye disorder or are sensitive to light.
- Do not use this product if you suffer from a condition that causes photosensitivity.
- Do not use this product for pregnant women, people with serious heart conditions, seriously ill patients, or other emergency patients. If in doubt, consult a physician first.
- Do not use this product near babies.
- Do not use the device for longer than 20 minutes per therapy session.
- If you lend the product to someone else, ensure that the product is provided together with the user manual.

Elektra

- Do not disassemble, repair, or modify the product. Doing so may cause damage. Warranty is void if you open or repair the product yourself. For safety, parts may only be replaced by authorized personnel. Always return the product to a service center for repairs to prevent damage or personal injury.
- Do not use the product near a bath, shower, swimming pool, or sink filled with water.
- Ensure the product does not come into contact with water. Do not immerse in water. If water does enter the product, switch it off immediately and contact Customer Service.
- Do not touch the product with wet hands.

Other

- Ensure that the buttons are not pressed unintentionally when plugging in or moving the product.
- Do not insert objects into the openings or slots of the product. Ensure that no foreign objects enter the product housing.

About the LED Face Mask

The LED Face Mask is a portable and rechargeable beauty device for home use, based on advanced LED phototherapy. The mask is made of high-quality medical-grade silicone, ensuring comfortable, safe, and easy use.

Thanks to its lightweight and foldable design, the mask is easy to carry. The adjustable nylon strap makes it suitable for different face shapes and head sizes.

The device is equipped with 3 LED chips per light point, ensuring maximum coverage and an even, effective treatment:


- 70 light points × 3 chips = 210 LEDs

This distribution ensures that all areas are illuminated evenly and intensively.

Each LED emits four different wavelengths: red light (630 nm), blue light (460 nm), yellow light (590 nm), and near-infrared light (850 nm).

The accompanying controller offers 7 adjustable color therapy modes, 5 intensity levels, and multiple adjustable treatment duration options. This allows the treatment to be easily tailored to personal preferences and skin needs.

Through LED phototherapy, the light penetrates deep into the skin, where it supports cellular activity and stimulates natural skin processes. This helps support collagen and connective tissue production, contributing to a firmer, more even, and fresher-looking skin. In addition, the device supports circulation, skin elasticity, and overall skin condition, which may contribute to a more even complexion, improved skin balance, and the softening of fine lines and pigmentation.

Receiving Shipment




- Check that the contents of the package match the packing slip when you receive the shipment. Notify your supplier if any parts are missing. If the equipment appears damaged, file a complaint with the carrier immediately and give the supplier a detailed description of the damage. Keep the damaged packaging to substantiate your complaint.

Package Contents

- Good4You LED Face Mask - Model MO8


Included Accessories

- 2× Adjustable fastening straps
- 2× Eye protectors
- 1× Charging cable
- 1× Remote control
- 1× User manual

Installation diagram

Attach the fastening straps
Connect the device to the
remote control


Use and Settings

- Clean your face thoroughly before using the LED Face Mask. Residues of make-up or skincare products may reduce the effectiveness of the light therapy.
- Fully charge the remote control before first use.
- Attach the LED Face Mask using the adjustable fastening straps.
 - » Thread the straps through the openings on both sides of the mask.
 - » Adjust the straps until the mask fits securely and comfortably on your face.
- Before first use, perform a light sensitivity test.
 - » Place the device on the inside of your arm and let it operate for 10 minutes.
 - » If redness, itching, or other skin reactions occur, stop using the device and consult a beautician before using it again.
 - » After use, the skin may appear temporarily red. This is a normal skin reaction and usually disappears within 1 to 2 hours.
- Always start the treatment at intensity level 1 or 2 so your skin can gradually adapt to the light therapy.
- Increase the intensity and treatment duration gradually, depending on your skin condition and personal comfort.
- The recommended usage time is 20 minutes per session.
- For optimal results, it is recommended to use the device 2 to 3 times per week.

Charging instructions

- ⓘ This product contains a built-in lithium-ion polymer battery. Do not use the device while charging, as this may cause damage to the battery and the device, which may result in improper operation
1. Fully charge the device before first use. To do so, connect the USB charging cable to the remote control and connect it to an external power source, such as a USB adapter or power bank.
 2. Indicator status: when the indicator light on the remote control is flashing, the device is charging. The charging time is approximately 4 to 5 hours. Once the indicator light turns off, charging is complete.
 3. When the battery is nearly empty, the device switches off automatically. Fully recharge the device before using it again

Explanation of Light Modes

- **Red light (630 nm)**

Supports natural cellular activity and skin renewal. May contribute to firmer skin, improved elasticity, and a smoother skin texture. Often used to help the skin appear more even, fresh, and youthful.

- **Blue light (460 nm)**

Helps cleanse the skin and may contribute to reducing impurities. Supports a calmer, clearer skin and is often used for unsettled or acne-prone skin.

- **Green light (520 nm)**

Has a balancing effect on the skin and may help regulate sebum production. Supports a more even skin tone and contributes to a calm and comfortable skin feel.

- **Purple light (630 nm + 460 nm – combination of red and blue)**

Combines the benefits of red and blue light. Supports skin recovery, helps calm the skin, and may contribute to a more even and healthy-looking complexion.

- **Yellow light (520 nm + 630 nm – combination of green and red)**

Supports skin balance and may help reduce visible skin irregularities. Often used to help the skin appear fresher and brighter.

- **Cyan light (520 nm + 460 nm – combination of green and blue)**

Supports skin metabolism and may contribute to a vital, fresh appearance. Often used to help energize tired-looking skin.

- **White light (460 nm + 520 nm + 630 nm – combination of blue, green, and red)**

Combines multiple wavelengths for broad skin support. Penetrates deeper into the skin and may help soften fine lines, improve skin texture, and reduce a tired appearance.

Technical data

Charging voltage	5V 0.5A
Battery capacity	1200mAh
Number of LEDs	210
Net weight	145g

Maintenance and Cleaning

Maintenance

- Turn off the product when not in use and unplug it.
- Keep the product dust-free.
- Make sure the product is not in contact with sharp objects that may damage the product.
- Do not fold the product sharply during storage, as this may cause damage to the internal wiring

Cleaning

Do not immerse the product in water or other liquids and keep it away from solvents and harsh cleaning agents.

- Turn off and unplug the product before cleaning. Clean the base and control panel of the product with a damp cloth dipped briefly in detergent and wipe the product dry with a dry, soft cloth. Do not use abrasives or thinners to clean the product.

Disposal and Recycling

- Think about the environment and contribute to a cleaner living environment! At the end of its service life, do not dispose of the product with normal household waste, but hand it in at a government designated collection point for recycling.

Service and Warranty

- If you need service or information regarding your product, please do not hesitate to contact us via email: info@good-4you.nl
- We provide one year warranty on our products. The customer must be able to show proof of purchase and the warranty must be validated within a period of one year from the date of purchase. Please note that the warranty does not apply: for damage caused by improper use or repair; for parts subject to wear and tear; for defects known to the customer at the time of purchase; for defects caused by the customer's own fault, for damage caused by third parties.

Extend warranty for FREE

1) Scan the QR code below to find Good4You's Facebook page and then 'Like' it.

2) Send a message with "Warranty" to us via Messenger or email info@good-4you.nl to extend the warranty from 1 year to 3 years.

Disclaimer

- Changes subject to change; specifications are subject to change without notice.

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© OnlineGoods