




Red Light Therapy Panel

Touch screen manual




START/PAUSE




TIME SETTING



Click to enter




Time setting, you can slide the color block to adjust it at will




You can also turn off the time setting

MODE SELECTION




Enter the lighting mode selection

Red 、 Nir or Blue




DIMMER

Click to enter brightness adjustment




Slide the color block to adjust the brightness, or click "+" or "-"




PULSE SETTING



Click on the numbers




Enter pulse mode



You can select or click to enter the desired pulse frequency



ADD AS DEFAULT MODE



Adding a custom mode


DEFAULT MODE



Click to enter the custom page

Slide the color block to the left to delete it


DEFAULT MODE




Click to enter the custom page

Long press the color block to enter the digital adjustment page


SYNCHRONOUS WORKING



Sync Press and hold for 5 seconds to enter the pairing interface



Click on the machine you want to pair with



Ok


SYNCHRONOUS WORKING



PRESET MODE



Muscle



Sleep



Mood

30:00

Muscle

Joint Pain

Sleep

Skin

Mood

Recovery



Joint Pain



Skin



Recovery

